



COMPANY PROFILE:

Company Name:

Sinus instruments .co

Establishment

1986

Business:

Manufacturers-cum-Exporters of
Surgical, Dental, Beauty Instruments
& all sort of Scissors

Telephones:

0092523301054

Cell:

+923049315806

Sinusinstruments@gmail.com



COMPANY PROFILE:

In 1986 a team of professional experts joined hands to form **Sinus Instruments .co** with the sole aim to Export High Quality Surgical and Dental Instruments. Under the banner of **Sinus Instruments** our products were advertised throughout the world by the help of Catalogs/ Brochures and by personal visits to different countries of the world. These efforts bore fruit and we managed to establish our credibility and began exporting quality surgical and dental products.

Sinus Instruments .co achieved tremendous success by virtue of constant dedicated efforts, devotion to duty and maintenance of quality control. This quality management system was properly implemented and is being successfully executed at all production stages. This is indeed a result of constant, professional and vigorous team effort undertaken by Buraq's most responsible, professional and dedicated staff members.

Keeping in view the universal quality standards, we are fortunate to obtain **ISO:13485:2012** Quality management System Certificate and also achieved the landmark certification of **ISO-9001: 2000**. Regular audit is carried out by experienced auditors.

It 's been almost 30 years of dedicated service which has earned trust and reputation to us in this field. In the meantime **Sinus** has recently ventured in a new field of Beauty Instruments with the brand name of **Sinus Instruments .co** We sincerely hope that this product will attract favorable attention from all over the world. Necessary information about this new product line is available in our beauty instruments catalogue, which can be provided on request.

It 's always been our endeavour to maintain a high standard of production, quality control and customer satisfaction. We finally offer our sincere regard and gratitude to all our valued clients for patronizing **Sinus Instruments .co** and helping us to achieve our basic target/goal.



Warranty

All instruments manufactured by **Sinus instruments .co** Will be replaced free of cost in case of any manufacturing/workmanship defect. However these instruments must be used as per the required International Standards/procedures

Quality Management System

All instruments are manufactured by **Sinus instruments .co** under strict quality Management / Regulations which conform to **ISO 9001:2000**

Special Note

This catalogue is sole property of **Sinus instruments .co** and any attempt to copy it by any means will be taken as a breach of the copyright Act and he or she shall be liable to prosecution under the Law.

This catalogue can be amended/modified by **Sinus instruments .co** without any prior notice.

Since all instruments in this catalogue are of International Pattern, **Sinus instruments .co** will not be responsible for any similarity and it will be viewed as coincidental.

We have tried to include maximum range of instruments in this catalogue and we guarantee to manufacture and supply any new item of the client's choice provided a sample/photograph along with necessary specification is given. We have also endeavoured clarity regarding sizes / patterns and shall accept any error on our part.

However as these instruments are handmade, marginal deviation in size and weight may occur which will not affect their performance.


Finally we suggest that these instruments must be handled only by professionals and **Sinus instruments .co** shall not be responsible for any mishap to human life or any physical organ by the misuse of these instruments by inexperienced surgeons.

BONE SCREWS

LOCKING & DCP



Technical Guidance



T.P.I = THREAD PER INCH
DIA = DIAMETER

DCP SCREWS

2.7 MM CORTICAL SCREW - SELF TAPPING



►REF. NO. SS	11-2708 to 11-2740		
►REF. NO. TT	11-2708T to 11-2740T		
►LENGTH IN MM	8 MM to 40 MM WITH 2 MM DIFF.		
►HEAD DIA:	5 MM	►HEXAGONAL SOCKET:	2 MM
►CORE DIA:	1.9 MM	►TAP DIA:	2.7 MM
►TPI:	1 MM	►DRILL BIT:	2 MM
►SHAFT DIA:	N/A	►GUIDE WIRE:	N/A
►THREAD DIA:	2.7 MM		

3.5 MM CORTICAL SCREW - 20TPI - SELF TAPPING



►REF. NO. SS	11-3510 to 11-3570		
►REF. NO. TT	11-3510T to 11-3570T		
►LENGTH IN MM	10 MM to 70 MM WITH 2 MM DIFF.		
►HEAD DIA:	6 MM	►HEXAGONAL SOCKET:	2.5 MM
►CORE DIA:	2.4 MM	►TAP DIA:	3.5 MM
►TPI:	1.25 MM	►DRILL BIT:	2.5 MM
►SHAFT DIA:	N/A	►GUIDE WIRE:	N/A
►THREAD DIA:	3.5 MM		

4.5 MM CORTICAL SCREW - 14TPI - SELF TAPPING



►REF. NO. SS	11-4510 to 11-4570		
►REF. NO. TT	11-4510T to 11-4570T		
►LENGTH IN MM	10 MM to 70 MM WITH 2 MM DIFF.		
►HEAD DIA:	8 MM	►HEXAGONAL SOCKET:	3.5 MM
►CORE DIA:	3 MM	►TAP DIA:	4.5 MM
►TPI:	1.75 MM	►DRILL BIT:	3.2 MM
►SHAFT DIA:	N/A	►GUIDE WIRE:	N/A
►THREAD DIA:	4.5 MM		

4.5 MM MALLEOLAR SCREW - SELF TAPPING



►REF. NO. SS	11-4525 to 11-4585		
►REF. NO. TT	N/A		
►LENGTH IN MM	25 MM to 85 MM WITH 5 MM DIFF.		
►HEAD DIA:	8 MM	►HEXAGONAL SOCKET:	3.5 MM
►CORE DIA:	3 MM	►TAP DIA:	N/A
►TPI:	1.75 MM	►GUIDE WIRE:	3.2 MM
►SHAFT DIA:	3 MM	►DRILL BIT:	N/A
►THREAD DIA:	4.5 MM		

DCP SCREWS

3.5 MM CANCELLOUS SCREW - PARTIAL THREADED



►REF. NO. SS	11-3510 to 11-3570		
►REF. NO. TT	N/A		
►LENGTH IN MM	10 MM to 30 MM WITH 2 MM DIFF. 30 MM to 70 MM WITH 5 MM DIFF.		
►HEAD DIA:	6 MM	►HEXAGONAL SOCKET:	2.5 MM
►CORE DIA:	1.9 MM	►TAP DIA:	N/A
►TPI:	1.75 MM	►DRILL BIT:	2 MM
►SHAFT DIA:	2.5 MM	►GUIDE WIRE:	N/A
►THREAD DIA:	3.5 MM		

3.5 MM CANCELLOUS SCREW - FULL THREADED



►REF. NO. SS	11-3510 to 11-3570		
►REF. NO. TT	N/A		
►LENGTH IN MM	10 MM to 30 MM WITH 2 MM DIFF.		
	30 MM to 70 MM WITH 5 MM DIFF.		
►HEAD DIA:	6 MM	►HEXAGONAL SOCKET:	2.5 MM
►CORE DIA:	1.9 MM	►TAP DIA:	3.5 MM
►TPI:	1.75 MM	►DRILL BIT:	2 MM
►SHAFT DIA:	N/A	►GUIDE WIRE:	N/A
►THREAD DIA:	3.5 MM		

4.0 MM CANCELLOUS SCREW - PARTIAL THREADED



►REF. NO. SS	11-4010 to 11-4070		
►REF. NO. TT	N/A		
►LENGTH IN MM	10 MM to 30 MM WITH 2 MM DIFF. 30 MM to 70 MM WITH 5 MM DIFF.		
►HEAD DIA:	6 MM	►HEXAGONAL SOCKET:	2.5 MM
►CORE DIA:	1.9 MM	►TAP DIA:	N/A
►TPI:	1.75 MM	►DRILL BIT:	2 MM
►SHAFT DIA:	2.5 MM	►GUIDE WIRE:	N/A
►THREAD DIA:	4 MM		

4.0 MM CANCELLOUS SCREW - FULL THREADED



►REF. NO. SS	11-4010 to 11-4070		
►REF. NO. TT	N/A		
►LENGTH IN MM	10 MM to 30 MM WITH 2 MM DIFF. 30 MM to 70 MM WITH 5 MM DIFF.		
►HEAD DIA:	6 MM	►HEXAGONAL SOCKET:	2.5 MM
►CORE DIA:	1.9 MM	►TAP DIA:	4 MM
►TPI:	1.75 MM	►DRILL BIT:	2 MM
►SHAFT DIA:	N/A	►GUIDE WIRE:	N/A
►THREAD DIA:	4 MM		

DCP SCREWS

6.5 MM CANCELLOUS SCREW - 16 MM THREADED

►REF. NO. SS	11-6525 to 11-65110		
►REF. NO. TT	N/A		
►LENGTH IN MM	25 MM to 110 MM WITH 5 MM DIFF.		
►HEAD DIA:	8 MM	►HEXAGONAL SOCKET:	3.5 MM
►CORE DIA:	3 MM	►TAP DIA:	N/A
►TPI:	2.8 MM	►DRILL BIT:	3.2 MM
►SHAFT DIA:	4.5 MM	►GUIDE WIRE:	N/A
►THREAD DIA:	6.5 MM		

6.5 MM CANCELLOUS SCREW - 32 MM THREADED

►REF. NO. SS	11-6540 to 11-65110		
►REF. NO. TT	N/A		
►LENGTH IN MM	40 MM to 110 MM WITH 5 MM DIFF.		
►HEAD DIA:	8 MM	►HEXAGONAL SOCKET:	3.5 MM
►CORE DIA:	3 MM	►TAP DIA:	N/A
►TPI:	2.8 MM	►DRILL BIT:	3.2 MM
►SHAFT DIA:	4.5 MM	►GUIDE WIRE:	N/A
►THREAD DIA:	6.5 MM		

6.5 MM CANCELLOUS SCREW - FULLY THREADED

►REF. NO. SS	11-6525 to 11-65110		
►REF. NO. TT	N/A		
►LENGTH IN MM	25 MM to 110 MM WITH 5 MM DIFF.		
►HEAD DIA:	8 MM	►HEXAGONAL SOCKET:	3.5 MM
►CORE DIA:	3 MM	►TAP DIA:	6.5 MM
►TPI:	2.8 MM	►DRILL BIT:	3.2 MM
►SHAFT DIA:	N/A	►GUIDE WIRE:	N/A
►THREAD DIA:	6.5 MM		

4.0 MM CANNULATED CANCELLOUS SCREW - PARTIALLY THREADED

►REF. NO. SS	11-4010 to 11-4070		
►REF. NO. TT	11-4010T to 11-4070T		
►LENGTH IN MM	10 MM to 70 MM WITH 5 MM DIFF.		
►HEAD DIA:	6 MM	►HEXAGONAL SOCKET:	2.5 MM
►CORE DIA:	2.6 MM	►TAP DIA:	N/A
►TPI:	1.75 MM	►GUIDE WIRE:	2.7 MM
►SHAFT DIA:	2.7 MM	►DRILL BIT:	1.2 MM
►THREAD DIA:	4 MM		

DCP SCREWS

4.0 MM CANNULATED CANCELLOUS SCREW - FULLY THREADED



▶REF. NO. SS	11-4010 to 11-4070		
▶REF. NO. TT	11-4010T to 11-4070T		
▶LENGTH IN MM	10 MM to 70 MM WITH 5 MM DIFF.		
▶HEAD DIA:	6 MM	▶HEXAGONAL SOCKET:	2.5 MM
▶CORE DIA:	2.6 MM	▶TAP DIA:	4 MM
▶TPI:	1.75 MM	▶DRILL BIT:	2.7 MM
▶SHAFT DIA:	N/A	▶GUIDE WIRE:	1.2 MM
▶THREAD DIA:	4 MM		

6.5 MM LARGE CANNULATED CANCELLOUS SCREW - 16 MM THREADED



▶REF. NO. SS	11-6525 to 11-65110		
▶REF. NO. TT	11-6525T to 11-65110T		
▶LENGTH IN MM	25 MM to 110 MM WITH 5 MM DIFF.		
▶HEAD DIA:	8 MM	▶HEXAGONAL SOCKET:	3.5 MM
▶CORE DIA:	4.4 MM	▶TAP DIA:	N/A
▶TPI:	2.8 MM	▶DRILL BIT:	4.5 MM
▶SHAFT DIA:	4.5 MM	▶GUIDE WIRE:	1.8 MM
▶THREAD DIA:	6.5 MM		

6.5 MM LARGE CANNULATED CANCELLOUS SCREW - 32 MM THREADED



▶REF. NO. SS	11-6540 to 11-65110		
▶REF. NO. TT	11-6540T to 11-65110T		
▶LENGTH IN MM	40 MM to 110 MM WITH 5 MM DIFF.		
▶HEAD DIA:	8 MM	▶HEXAGONAL SOCKET:	3.5 MM
▶CORE DIA:	4.4 MM	▶TAP DIA:	N/A
▶TPI:	2.8 MM	▶DRILL BIT:	4.5 MM
▶SHAFT DIA:	4.5 MM	▶GUIDE WIRE:	1.8 MM
▶THREAD DIA:	6.5 MM		

6.5 MM LARGE CANNULATED CANCELLOUS SCREW - FULLY THREADED



▶REF. NO. SS	11-6530 to 11-65115		
▶REF. NO. TT	11-6530T to 11-65115T		
▶LENGTH IN MM	30 MM to 115 MM WITH 5 DIFF.		
▶HEAD DIA:	8 MM	▶HEXAGONAL SOCKET:	3.5 MM
▶CORE DIA:	4.4 MM	▶TAP DIA:	6.5 MM
▶TPI:	2.8 MM	▶DRILL BIT:	4.5 MM
▶SHAFT DIA:	N/A	▶GUIDE WIRE:	1.8 MM
▶THREAD DIA:	6.5 MM		

DCP & LCP SCREWS

12.5 MM DHS SCREW WITH COMPRESSION SCREW



▶REF. NO. SS	11-1250 to 1123-SS-12110		
▶REF. NO. TT	N/A		
▶LENGTH IN MM	50 MM to 110 MM WITH 5 MM DIFF.		
▶HEAD DIA:	N/A	▶HEXAGONAL SOCKET:	3.5 MM
▶CORE DIA:	8 MM	▶TAP DIA:	12.7 MM
▶TPI:	3.2 MM	▶DRILL BIT:	8 MM
▶SHAFT DIA:	8 MM	▶GUIDE WIRE:	2.5 MM
▶THREAD DIA:	12.7 MM		

3.5 MM HERBERT SCREW (HEADLESS)



▶REF. NO. SS	N/A		
▶REF. NO. TT	11-3510T to 11-3530T		
▶LENGTH IN MM	10 MM to 30 MM WITH 2 MM DIFF.		
▶HEAD DIA:	4.5 MM	▶HEXAGONAL SOCKET:	2.0 MM
▶CORE DIA:	2.7 MM	▶TAP DIA:	3.5 MM
▶TPI:	1.25 MM	▶DRILL BIT:	2.7 MM
▶SHAFT DIA:	2.7 MM	▶GUIDE WIRE:	1.2 MM
▶THREAD DIA:	3.5 MM		

2.4 MM LCP SCREW



▶REF. NO. SS	14-2408 to 14-2440		
▶REF. NO. TT	14-2408T to 14-2440T		
▶LENGTH IN MM	8 MM to 40 MM WITH 2 MM DIFF.		
▶HEAD DIA:	3.4 MM	▶HEXAGONAL SOCKET:	2 MM
▶CORE DIA:	1.9 MM	▶TAP DIA:	N/A
▶TPI:	0.6 MM	▶DRILL BIT:	1.8 MM
▶SHAFT DIA:	N/A	▶GUIDE WIRE:	N/A
▶THREAD DIA:	2.4 MM		

2.7 MM LCP SCREW



▶REF. NO. SS	14-2708 to 14-2740		
▶REF. NO. TT	14-2708T to 14-2740T		
▶LENGTH IN MM	8 MM to 40 MM WITH 2 MM DIFF.		
▶HEAD DIA:	3.4 MM	▶HEXAGONAL SOCKET:	2 MM
▶CORE DIA:	2.1 MM	▶TAP DIA:	N/A
▶TPI:	0.6 MM	▶DRILL BIT:	2 MM
▶SHAFT DIA:	N/A	▶GUIDE WIRE:	N/A
▶THREAD DIA:	2.7 MM		

LCP SCREWS

3.5 MM LCP SCREW



►REF. NO. SS	14-3510 to 14-3570		
►REF. NO. TT	14-3510T to 14-3570T		
►LENGTH IN MM	10 MM to 60 MM WITH 2 MM DIFF. 60 MM to 70 MM WITH 5 MM DIFF.		
►HEAD DIA:	5 MM	►HEXAGONAL SOCKET:	2.5 MM
►CORE DIA:	2.9 MM	►TAP DIA:	N/A
►TPI:	0.8 MM	►DRILL BIT:	2.8 MM
►SHAFT DIA:	N/A	►GUIDE WIRE:	N/A
►THREAD DIA:	3.5 MM		

5.0 MM LCP SCREW



►REF. NO. SS	14-5014 to 14-5090		
►REF. NO. TT	14-5014T to 14-5090T		
►LENGTH IN MM	14 MM to 60 MM WITH 2 MM DIFF.		
	60 MM to 90 MM WITH 5 MM DIFF.		
►HEAD DIA:	6.5 MM	►HEXAGONAL SOCKET:	3.5 MM
►CORE DIA:	4.4 MM	►TAP DIA:	N/A
►TPI:	1 MM	►DRILL BIT:	4.3 MM
►SHAFT DIA:	N/A	►GUIDE WIRE:	N/A
►THREAD DIA:	5 MM		

4.0 MM LCP CANCELLOUS SCREW PARTIALLY THREADED



►REF. NO. SS	14-4010 to 14-4070		
►REF. NO. TT	14-4010T to 14-4070T		
►LENGTH IN MM	10 MM to 30 MM WITH 2 MM DIFF. 30 MM to 70 MM WITH 5 MM DIFF.		
►HEAD DIA:	5 MM	►HEXAGONAL SOCKET:	2.5 MM
►CORE DIA:	2.4 MM	►TAP DIA:	N/A
►TPI:	1.75 MM	►DRILL BIT:	2.5 MM
►SHAFT DIA:	2.5 MM	►GUIDE WIRE:	N/A
►THREAD DIA:	4 MM		

4.0 MM LCP CANCELLOUS SCREW FULLY THREADED



►REF. NO. SS	14-4010 to 14-4070		
►REF. NO. TT	14-4010T to 14-4070T		
►LENGTH IN MM	10 MM to 30 MM WITH 2 MM DIFF. 30 MM to 70 MM WITH 5 MM DIFF.		
►HEAD DIA:	5 MM	►HEXAGONAL SOCKET:	2.5 MM
►CORE DIA:	2.4 MM	►TAP DIA:	N/A
►TPI:	1.75 MM	►DRILL BIT:	2.5 MM
►SHAFT DIA:	N/A	►GUIDE WIRE:	N/A
►THREAD DIA:	4 MM		

LCP SCREWS

5.0 MM LCP CANCELLOUS SCREW 32 MM THREADED



▶REF. NO. SS	14-5040 to 14-5090		
▶REF. NO. TT	14-5040T to 14-5090T		
▶LENGTH IN MM	40 MM to 90 MM WITH 5 MM DIFF.		
▶HEAD DIA:	6.5 MM	▶HEXAGONAL SOCKET:	3.5 MM
▶CORE DIA:	3.4 MM	▶TAP DIA:	N/A
▶TPI:	2.8 MM	▶DRILL BIT:	3.5 MM
▶SHAFT DIA:	4.5 MM	▶GUIDE WIRE:	N/A
▶THREAD DIA:	5 MM		

5.0 MM LCP CANCELLOUS SCREW FULLY THREADED



▶REF. NO. SS	14-5025 to 14-5090		
▶REF. NO. TT	14-5025T to 14-5090T		
▶LENGTH IN MM	25 MM to 90 MM WITH 5 MM DIFF.		
▶HEAD DIA:	6.5 MM	▶HEXAGONAL SOCKET:	3.5 MM
▶CORE DIA:	3.4 MM	▶TAP DIA:	N/A
▶TPI:	2.8 MM	▶DRILL BIT:	3.5 MM
▶SHAFT DIA:	N/A	▶GUIDE WIRE:	N/A
▶THREAD DIA:	5 MM		

6.5 MM LCP CANCELLOUS SCREW 16 MM THREADED



▶REF. NO. SS	14-6525 to 14-65100		
▶REF. NO. TT	N/A		
▶LENGTH IN MM	25 MM to 100 MM WITH 5 DIFF.		
▶HEAD DIA:	8 MM	▶HEXAGONAL SOCKET:	3.5 MM
▶CORE DIA:	3.4 MM	▶TAP DIA:	N/A
▶TPI:	2.8 MM	▶DRILL BIT:	3.5 MM
▶SHAFT DIA:	4.5 MM	▶GUIDE WIRE:	N/A
▶THREAD DIA:	6.5 MM		

6.5 MM LCP CANCELLOUS SCREW 32 MM THREADED



▶REF. NO. SS	14-6540 to 14-65100		
▶REF. NO. TT	N/A		
▶LENGTH IN MM	40 MM to 100 MM WITH 5 MM DIFF.		
▶HEAD DIA:	8 MM	▶HEXAGONAL SOCKET:	3.5 MM
▶CORE DIA:	3.4 MM	▶TAP DIA:	N/A
▶TPI:	2.8 MM	▶GUIDE WIRE:	3.5 MM
▶SHAFT DIA:	4.5 MM	▶DRILL BIT:	N/A
▶THREAD DIA:	6.5 MM		

LCP SCREWS

6.5 MM LCP CANCELLOUS SCREW FULLY THREADED



▶REF. NO. SS	14-6525 to 14-6590		
▶REF. NO. TT	N/A		
▶LENGTH IN MM	25 MM to 90 MM WITH 5 MM DIFF.		
▶HEAD DIA:	8 MM	▶HEXAGONAL SOCKET:	3.5 MM
▶CORE DIA:	3.4 MM	▶TAP DIA:	N/A
▶TPI:	2.8 MM	▶DRILL BIT:	3.5 MM
▶SHAFT DIA:	N/A	▶GUIDE WIRE:	N/A
▶THREAD DIA:	6.5 MM		

3.5 MM LCP SPACER



▶REF. NO. SS	14-0035		
▶REF. NO. TT	14-0050T		
▶LENGTH IN MM	Ø 3.5 MM X 2.0 MM		
▶HEAD DIA:	5.0 MM	▶HEXAGONAL SOCKET:	2.5 MM
▶CORE DIA:	N/A	▶TAP DIA:	N/A
▶TPI:	N/A	▶DRILL BIT:	2.8 MM
▶SHAFT DIA:	3.5 MM	▶GUIDE WIRE:	N/A
▶THREAD DIA:	N/A		

5.0 MM LCP SPACER



▶REF. NO. SS	14-0050		
▶REF. NO. TT	14-0050T		
▶LENGTH IN MM	Ø 5.0 MM X 2.0 MM		
▶HEAD DIA:	6.5 MM	▶HEXAGONAL SOCKET:	3.5 MM
▶CORE DIA:	N/A	▶TAP DIA:	N/A
▶TPI:	N/A	▶DRILL BIT:	4.3 MM
▶SHAFT DIA:	5 MM	▶GUIDE WIRE:	N/A
▶THREAD DIA:	N/A		

13.0 / 6.5 MM SPIKED WASHER FOR 6.5 MM CANCELLOUS SCREW.



▶REF. NO. SS	11-0027		
▶REF. NO. TT	11-0027T		
▶INNER DIA:	6.6 MM	▶OUTER DIA:	13 MM
▶THICKNESS	5 MM		

LCP SCREWS

13.0 / 6.5 MM WASHER FOR 6.5 MM CANCELLOUS SCREW



►REF. NO. SS 11-0035
►REF. NO. TT 11-0035T

►INNER DIA: 6.6 MM ►OUTER DIA: 13 MM
►THICKNESS 1 MM

7.0 / 4.0 MM WASHER FOR 4.0 MM CANCELLOUS SCREW.



►REF. NO. SS 11-0040
►REF. NO. TT 11-0040T

►INNER DIA: 4.2 MM ►OUTER DIA: 7 MM
►THICKNESS 1 MM

7.0 / 4.0 MM SPIKED WASHER FOR 4.0 MM CANCELLOUS SCREW.



►REF. NO. SS 11-0045
►REF. NO. TT 11-0045T

►INNER DIA: 4.2 MM ►OUTER DIA: 7 MM
►THICKNESS 5 MM

10.0 / 4.5 MM WAHSER FOR MALLEOR SCREW.



►REF. NO. SS 11-0065
►REF. NO. TT 11-0065T

►INNER DIA: 4.6 MM ►OUTER DIA: 10 MM
►THICKNESS 1 MM

DCP PLATES

DCP PLATES



DCP PLATES

2.7 MM DCP PLATE



►REF. NO. SS	12-2702 to 12-2712
►REF. NO. TT	N/A
►HOLES	2, 4, 5, 6, 7, 8, 9, 10, 11, 12

3.5 MM DCP SMALL



►REF. NO. SS	12-3504 to 12-3516
►REF. NO. TT	N/A
►HOLES	4, 5, 6, 7, 8, 9, 10, 12, 14, 16

4.5 MM DCP NARROW



►REF. NO. SS	12-4504 to 12-4518
►REF. NO. TT	N/A
►HOLES	4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 18

4.5 MM DCP BROAD



►REF. NO. SS	12-4506 to 12-4520
►REF. NO. TT	N/A
►HOLES	6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 18, 20

DCP PLATES

3.5 MM LC-DCP SMALL



►REF. NO. SS	12-3504 to 12-3516
►REF. NO. TT	N/A
►HOLES	4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16

4.5 MM LC-DCP NARROW



►REF. NO. SS	12-4504 to 12-4518
►REF. NO. TT	N/A
►HOLES	4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 18

4.5 MM LC-DCP BROAD



►REF. NO. SS	12-4505 to 12-4516
►REF. NO. TT	N/A
►HOLES	5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16

3.5 MM 1/3 TUBULAR PLATE



►REF. NO. SS	12-3503 to 12-3512
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

DCP PLATES

4.5 MM SEMI TUBULAR PLATE



►REF. NO. SS	12-4503 to 12-4512
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

3.5 MM STRAIGHT RECONSTRUCTION PLATE - DYNAMIC HOLES



►REF. NO. SS	12-3505 to 12-3516
►REF. NO. TT	N/A
►HOLES	5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 16

4.5 MM STRAIGHT RECONSTRUCTION PLATE - DYNAMIC HOLES



►REF. NO. SS	12-4505 to 12-4514
►REF. NO. TT	N/A
►HOLES	5, 6, 7, 8, 9, 10, 11, 12, 13, 14

3.5 MM STRAIGHT RECONSTRUCTION PLATE - ROUND HOLES (THICKNESS 1.6 MM)



►REF. NO. SS	12-3505 to 12-3524
►REF. NO. TT	N/A
►HOLES	5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24

DCP PLATES

4.5 MM STRAIGHT RECONSTRUCTION PLATE – ROUND HOLES (THICKNESS 1.6 MM)



►REF. NO. SS	12-4505 to 12-4524
►REF. NO. TT	N/A
►HOLES	5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24

3.5 MM SMALL T-PLATE, RIGHT ANGLE.



►REF. NO. SS	12-3503 to 12-3508
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8

3.5 MM SMALL T-PLATE, OBLIQUE ANGLE



►REF. NO. SS	12-3503 to 12-3508
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8

4.5 MM T-PLATE



►REF. NO. SS	12-4504 to 12-4512
►REF. NO. TT	N/A
►HOLES	4, 5, 6, 7, 8, 9, 10, 11, 12

DCP PLATES

4.5 MM T-BUTTRESS PLATE



►REF. NO. SS	12-4504 to 12-4512
►REF. NO. TT	N/A
►HOLES	4, 5, 6, 7, 8, 9, 10, 11, 12

4.5 MM L-BUTTRESS PLATE



RIGHT SIDE

►REF. NO. SS	12-4504 to 12-4512
►REF. NO. TT	N/A

LEFT SIDE

►REF. NO. SS	12-4504 to 12-4512
►REF. NO. TT	N/A
►HOLES	4, 5, 6, 7, 8, 9, 10, 11, 12

4.5 MM LATERAL TIBIAL BUTTRESS PLATE DCP HOLES (HOCKEY)



RIGHT SIDE

►REF. NO. SS	12-4505 to 12-4513
►REF. NO. TT	N/A

LEFT SIDE

►REF. NO. SS	12-4505 to 12-4513
►REF. NO. TT	N/A
►HOLES	5, 7, 9, 11, 13

3.5MM CLOVERLEAF PLATE



►REF. NO. SS	12-3503 to 12-3510
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10


DCP PLATES

3.5 MM SPOON PLATE



►REF. NO. SS	12-3504 to 12-3506
►REF. NO. TT	N/A
►HOLES	4, 5, 6

3.5 MM Y RECONSTRUCTION PLATE



►REF. NO. SS	12-3503 to 12-3507
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7

125° DHS PLATE LONG BARREL DYNAMIC HOLES



►REF. NO. SS	12-2502 to 12-2514
►REF. NO. TT	N/A
►HOLES	2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14

130° DHS PLATE LONG BARREL DYNAMIC HOLES



►REF. NO. SS	12-3002 to 12-3014
►REF. NO. TT	N/A
►HOLES	2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14

DCP PLATES

135° DHS PLATE LONG BARREL DYNAMIC HOLES



►REF. NO. SS	12-3502 to 12-3520
►REF. NO. TT	N/A
►HOLES	2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18, 20

140° DHS PLATE LONG BARREL DYNAMIC HOLES



►REF. NO. SS	12-4002 to 12-4016
►REF. NO. TT	N/A
►HOLES	2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16

145° DHS PLATE LONG BARREL DYNAMIC HOLES



►REF. NO. SS	12-4502 to 12-4516
►REF. NO. TT	N/A
►HOLES	2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16

150° DHS PLATE LONG BARREL DYNAMIC HOLES



►REF. NO. SS	12-5002 to 12-5016
►REF. NO. TT	N/A
►HOLES	2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16

DCP PLATES

125° DHS PLATE SHORT BARREL DYNAMIC HOLES



►REF. NO. SS	12-2503 to 12-2512
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

130° DHS PLATE SHORT BARREL DYNAMIC HOLES



►REF. NO. SS	12-3003 to 12-3012
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

135° DHS PLATE SHORT BARREL DYNAMIC HOLES



►REF. NO. SS	12-3503 to 12-3512
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

140° DHS PLATE SHORT BARREL DYNAMIC HOLES



►REF. NO. SS	12-4003 to 12-4012
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

DCP PLATES

145° DHS PLATE SHORT BARREL DYNAMIC HOLES



►REF. NO. SS	12-4503 to 12-4512
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

150° DHS PLATE SHORT BARREL DYNAMIC HOLES



►REF. NO. SS	12-5003 to 12-5012
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

125° DHS PLATE LONG BARREL ROUND HOLES



►REF. NO. SS	12-2502 to 12-2514
►REF. NO. TT	N/A
►HOLES	2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14

130° DHS PLATE LONG BARREL ROUND HOLES



►REF. NO. SS	12-3002 to 12-3014
►REF. NO. TT	N/A
►HOLES	2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14

DCP PLATES

135° DHS PLATE LONG BARREL ROUND HOLES



►REF. NO. SS	12-3502 to 12-3520
►REF. NO. TT	N/A
►HOLES	2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18, 20

140° DHS PLATE LONG BARREL ROUND HOLES



►REF. NO. SS	12-4002 to 12-4016
►REF. NO. TT	N/A
►HOLES	2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16

145° DHS PLATE LONG BARREL ROUND HOLES



►REF. NO. SS	12-4502 to 12-4516
►REF. NO. TT	N/A
►HOLES	2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16

150° DHS PLATE LONG BARREL ROUND HOLES



►REF. NO. SS	12-5004 to 12-5016
►REF. NO. TT	N/A
►HOLES	2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16

DCP PLATES

125° DHS PLATE SHORT BARREL ROUND HOLES



►REF. NO. SS	12-2503 to 12-2512
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

130° DHS PLATE SHORT BARREL ROUND HOLES



►REF. NO. SS	12-3003 to 12-3012
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

135° DHS PLATE SHORT BARREL ROUND HOLES



►REF. NO. SS	12-3503 to 12-3512
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

140° DHS PLATE SHORT BARREL ROUND HOLES



►REF. NO. SS	12-4003 to 12-4012
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

DCP PLATES

145° DHS PLATE SHORT BARREL ROUND HOLES



►REF. NO. SS	12-4503 to 12-4512
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

150° DHS PLATE SHORT BARREL ROUND HOLES



►REF. NO. SS	12-5003 to 12-5012
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

95° DCS PLATE DYNAMIC HOLES



►REF. NO. SS	12-9504 to 12-9518
►REF. NO. TT	N/A
►HOLES	4, 5, 6, 7, 8, 9, 10, 11, 12, 14, 16, 18

95° DCS PLATE ROUND HOLES



►REF. NO. SS	12-9504 to 12-9518
►REF. NO. TT	N/A
►HOLES	4, 5, 6, 7, 8, 9, 10, 11, 12, 14, 16, 18

DCP PLATES

95° DCP CONDYLAR BLADE PLATE



►REF. NO. SS	12-9505 to 12-9512
►REF. NO. TT	N/A
►HOLES	5, 7, 9, 12 - 50 MM, 60 MM, 70 MM, 80 MM

130° LCP DHS PLATE LONG BARREL DC HOLES



►REF. NO. SS	12-3003 to 12-3012
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

135° LCP DHS PLATE LONG BARREL DC HOLES



►REF. NO. SS	12-3503 to 12-3512
►REF. NO. TT	N/A
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12


INTRAMEDULLARY NAIL



INTRAMEDULLARY NAIL

Technical

Guidance



INTRAMEDULLARY NAIL

HUMERUS NAIL – CANNULATED

►DIA Ø	►REF. NO. SS
6.5 MM	16-6518 to 16-6530
7.0 MM	16-7018 to 16-7030
7.5 MM	16-7518 to 16-7530
►LENGTH IN MM 180, 200, 220, 240, 260, 280, 300	

TIBIA NAIL – CANNULATED

►DIA Ø	►REF. NO. SS	►REF. NO. TT
8.0 MM	16-8028 to 16-8038	16-8028T to 16-8038 T
9.0 MM	16-9028 to 16-9038	16-9028T to 16-9038 T
10 MM	16-1028 to 16-1038	16-1028T to 16-1038 T
11 MM	16-1128 to 16-1138	16-1128T to 16-1138 T
►LENGTH IN MM 280, 300, 320, 340, 360, 380		

EXPERT TIBIA NAIL – CANNULATED

►DIA Ø	►REF. NO. SS
8.0 MM	16-8028 to 16-8038
9.0 MM	16-9028 to 16-9038
10 MM	16-1028 to 16-1038
11 MM	16-1128 to 16-1138
►LENGTH IN MM 280, 300, 320, 340, 360, 380	

FEMUR NAIL – CANNULATED

►DIA Ø	►REF. NO. SS	►REF. NO. TT
9.0 MM	16-9032 to 16-9044	16-9032T to 16-9044T
10 MM	16-1032 to 16-1044	16-1032T to 16-1044T
11 MM	16-1132 to 16-1144	16-1132T to 16-1144T
12 MM	16-1232 to 16-1244	16-1232T to 16-1244T
►LENGTH IN MM 320, 340, 360, 380, 400, 420, 440		

INTRAMEDULLARY NAIL

130° SHORT PROXIMAL FEMORAL NAIL

►DIA Ø	►REF. NO. SS	►REF. NO. TT
9.0 MM	16-9025	16-9025T
10 MM	16-1025	16-1025T
11 MM	16-1125	16-1125T
12 MM	16-1225	16-1225T
►LENGTH IN MM 250 MM		

135° SHORT PROXIMAL FEMORAL NAIL

►DIA Ø	►REF. NO. SS	►REF. NO. TT
9.0 MM	16-9025	16-9025T
10 MM	16-1025	16-1025T
11 MM	16-1125	16-1125T
12 MM	16-1225	16-1225T
►LENGTH IN MM 250 MM		

130° TROCHANTERIC FEMUR NAIL

►DIA Ø	►REF. NO. SS	►REF. NO. TT
9.0 MM	16-9018	16-9018T
10 MM	16-1018	16-1018T
11 MM	16-1118	16-1118T
12 MM	16-1218	16-1218T
►LENGTH IN MM 250 MM		

135° TROCHANTERIC FEMUR NAIL

►DIA Ø	►REF. NO. SS	►REF. NO. TT
9.0 MM	16-9018	16-9018T
10 MM	16-1018	16-1018T
11 MM	16-1118	16-1118T
12 MM	16-1218	16-1218T
►LENGTH IN MM 180 MM		

INTRAMEDULLARY NAIL

130° - LONG PROXIMAL FEMORAL NAIL

RIGHT SIDE

►DIA Ø	►REF. NO. SS	►REF. NO. TT
9.0 MM	16-9034 to 16-9044	16-9034 T to 16-9044 T
10 MM	16-1034 to 16-1044	16-1034 T to 16-1044 T
11 MM	16-1134 to 16-1144	16-1134 T to 16-1144 T
►LENGTH IN MM 340, 360, 380, 400, 420, 440		



130° - LONG PROXIMAL FEMORAL NAIL

LEFT SIDE

►DIA Ø	►REF. NO. SS	►REF. NO. TT
9.0 MM	16-9034 to 16-9044	16-9034T to 16-9044T
10 MM	16-1034 to 16-1044	16-1034T to 16-1044T
11 MM	16-1134 to 16-1144	16-1134T to 16-1144T
►LENGTH IN MM 340, 360, 380, 400, 420, 440		



135° - LONG PROXIMAL FEMORAL NAIL

RIGHT SIDE

►DIA Ø	►REF. NO. SS	►REF. NO. TT
9.0 MM	16-9034 to 16-9044	16-9034T to 16-9044T
10 MM	16-1034 to 16-1044	16-1034T to 16-1044T
11 MM	16-1134 to 16-1144	16-1134T to 16-1144T
12 MM	16-1234 to 16-1244	16-1234T to 16-1244T
►LENGTH IN MM 340, 360, 380, 400, 420, 440		



135° - LONG PROXIMAL FEMORAL NAIL

LEFT SIDE

►DIA Ø	►REF. NO. SS	►REF. NO. TT
9.0 MM	16-9034 to 16-9044	16-9034 T to 16-9044 T
10 MM	16-1034 to 16-1044	16-1034 T to 16-1044 T
11 MM	16-1134 to 16-1144	16-1134 T to 16-1144 T
12 MM	16-1234 to 16-1244	16-1234 T to 16-1244 T
►LENGTH IN MM 340, 360, 380, 400, 420, 440		



INTRAMEDULLARY NAIL

SUPRACODYLAR DISTAL FEMUR NAIL



►DIA Ø	►REF. NO. SS
9.0 MM	16-9015 to 16-9030
10 MM	16-1015 to 16-1030
11 MM	16-1115 to 16-1130
12 MM	16-1215 to 16-1230
►LENGTH IN MM 150, 200, 250, 300	

ENDER'S NAIL



►DIA Ø	►REF. NO. SS	►LENGTH IN MM
2.0 MM	16-2020 to 16-2030	200 to 300 MM WITH 10 MM DIFF.
2.5 MM	16-2520 to 16-2530	200 to 300 MM WITH 10 MM DIFF.
3.0 MM	16-3021 to 16-3030	210 to 300 MM WITH 10 MM DIFF.
3.5 MM	16-3521 to 16-3530	210 to 300 MM WITH 10 MM DIFF.
4.0 MM	16-4021 to 16-4040	210 to 400 MM WITH 10 MM DIFF.
4.5 MM	16-4521 to 16-4540	210 to 400 MM WITH 10 MM DIFF.

TITANIUM ELASTIC NAIL - TENS NAIL



►DIA Ø	►REF. NO. TT	►LENGTH IN MM
2.0 MM	16-2044T	440
2.5 MM	16-2544T	440
3.0 MM	16-3044T	440
3.5 MM	16-3544T	440
4.0 MM	16-4044T	440
4.5 MM	16-4544T	440

K - NAIL



►DIA Ø	►REF. NO. SS	►LENGTH IN MM
7.0 MM	16-7032 to 16-7044	320, 340, 360, 380, 400, 420, 440
8.0 MM	16-8032 to 16-8044	320, 340, 360, 380, 400, 420, 440
9.0 MM	16-9032 to 16-9044	320, 340, 360, 380, 400, 420, 440
10.0 MM	16-1032 to 16-1044	320, 340, 360, 380, 400, 420, 440
11.0 MM	16-1132 to 16-1144	320, 340, 360, 380, 400, 420, 440
12.0 MM	16-1232 to 16-1244	320, 340, 360, 380, 400, 420, 440

INTRAMEDULLARY NAIL

END CAP - TIBIA/FEMUR/PFN/SUPRA



►DIA Ø	►REF. NO. SS	►REF. NO. TT
0 MM	16-0000	16-0000T
5 MM	16-0500	16-0500T
10 MM	16-1000	16-1000T
15 MM	16-1500	16-1500T
►LENGTH IN MM 0, 5, 10, 15		

INTERLOCKING BOLT, FOR HUMERUS NAIL



►DIA Ø	►REF. NO. SS	►REF. NO. TT
3.0 MM	17-0014 to 17-0040	
►LENGTH IN MM 14 MM to 40 MM WITH 2 MM DIFF.		

INTERLOCKING BOLT



►DIA Ø	►REF. NO. SS	►REF. NO. TT
3.9 MM	17-4024 to 17-4060	17-4024T to 17-4060T
►LENGTH IN MM 24 MM to 60 MM WITH 2 MM DIFF.		
4.9 MM	17-5024 to 17-5080	17-5024T to 17-5080T
►LENGTH IN MM 24 MM to 80 MM WITH 2 MM DIFF.		

6.4 MM HIP PIN (FOR PFN NAIL)



►DIA Ø	►REF. NO. SS	►REF. NO. TT
6.4 MM	17-6440 to 17-64110	17-6440T to 17-64110T
►LENGTH IN MM 40 MM to 110 MM WITH 5 MM DIFF.		

INTRAMEDULLARY NAIL

8.0 MM FEMORAL NECK SCREW (FOR PFN NAIL)

►DIA Ø	►REF. NO. SS	►REF. NO. TT
8.0 MM	17-8040 to 17-80110	17-8040T to 17-80110T
►LENGTH IN MM 40 MM to 110 MM WITH 5 MM DIFF.		

6.4 MM SOLID INTERLOCKING BOLT (FOR SUPRACONDYLAR)

►DIA Ø	►REF. NO. SS
6.4 MM	17-6440 to 17-64110
►LENGTH IN MM 40 MM to 110 MM WITH 5 MM DIFF.	

4.0 MM INTERLOCKING BOLT PARTIAL THREAD (FOR HUMERUS NAIL)

►DIA Ø	►REF. NO. SS
4.0 MM	17-4020 to 17-4050
►LENGTH IN MM 20 MM to 50 MM WITH 5 MM DIFF.	

SQUARE NAIL FOR ULNA

►DIA Ø	►REF. NO. SS	►LENGTH IN MM
2.0 MM	16-2020 to 16-2030	200, 220, 240, 260, 280, 300
2.5 MM	16-2520 to 16-2530	200, 220, 240, 260, 280, 300
3.0 MM	16-3020 to 16-3030	200, 220, 240, 260, 280, 300
3.5 MM	16-3520 to 16-3530	200, 220, 240, 260, 280, 300

INTRAMEDULLARY NAIL

SQUARE NAIL FOR RADIUS

►DIA Ø	►REF. NO. SS	►LENGTH IN MM
2.0 MM	16-2020 to 16-2030	200, 220, 240, 260, 280, 300
2.5 MM	16-2520 to 16-2530	200, 220, 240, 260, 280, 300
3.0 MM	16-3020 to 16-3030	200, 220, 240, 260, 280, 300
3.5 MM	16-3520 to 16-3530	200, 220, 240, 260, 280, 300


RUSH NAIL

►DIA Ø	►REF. NO. SS	►LENGTH IN MM
2.0 MM	16-2020 to 16-2030	200, 220, 240, 260, 280, 300
2.5 MM	16-2520 to 16-2530	200, 220, 240, 260, 280, 300
3.0 MM	16-3020 to 16-3030	200, 220, 240, 260, 280, 300
3.5 MM	16-3520 to 16-3530	200, 220, 240, 260, 280, 300
4.0 MM	16-4020 to 16-4030	200, 220, 240, 260, 280, 300

JOINT PROSTHESIS IMPLANTS



PROSTHESIS



PROSTHESIS

BIPOLAR PROSTHESIS - FENESTRATED - NEW TYPE - STERILISED

- ▶REF. NO. SS 26-3700 to 26-5500
- ▶REF. NO. TT N/A
- ▶LENGTH IN MM 37, 39, 41, 43, 45, 47, 51, 53, 55

BIPOLAR PROSTHESIS - NON FENESTRATED - NEW TYPE - STERILISED

- ▶REF. NO. SS 26-3700 to 26-5500
- ▶REF. NO. TT N/A
- ▶LENGTH IN MM 37, 39, 41, 43, 45, 47, 51, 53, 55

A.M. HIP PROSTHESIS - STERILISED

- ▶REF. NO. SS 26-3700 to 26-5500
- ▶REF. NO. TT N/A
- ▶LENGTH IN MM 37, 38, 39, 40, 41, 42, 43, 44, 45, 46, 47
48, 49, 50, 51, 52, 53, 54, 55

A.M. HIP PROSTHESIS - NON- STERILISED

- ▶REF. NO. SS 26-3700 to 26-5500
- ▶REF. NO. TT N/A
- ▶LENGTH IN MM 37, 38, 39, 40, 41, 42, 43, 44, 45, 46, 47
48, 49, 50, 51, 52, 53, 54, 55



PROSTHESIS

MODULAR BIPOLAR CUP

►REF. NO. SS	26-3700 to 26-5500
►REF. NO. TT	N/A
►LENGTH IN MM	
SMALL	37, 39, 41, 43, 45, 47, 49, 51, 53, 54, 55
MEDIUM	37, 39, 41, 43, 45, 47, 49, 51, 53, 54, 55
LARGE	37, 39, 41, 43, 45, 47, 49, 51, 53, 54, 55



MODULAR BIPOLAR STEM

►REF. NO. SS	26-0000 to 26-0000
►REF. NO. TT	N/A
►LENGTH	SMALL, MEDIUM & LARGE



THOMPSON HIP PROSTHESIS - STERILISED

►REF. NO. SS	26-0037 to 26-0055
►REF. NO. TT	N/A
►LENGTH IN MM	37, 38, 39, 40, 41, 42, 43, 44, 45, 46, 47 48, 49, 50, 51, 52, 53, 54, 55



THOMPSON HIP PROSTHESIS NON- STERILISED

►REF. NO. SS	26-0037 to 26-0055
►REF. NO. TT	N/A
►LENGTH IN MM	37, 38, 39, 40, 41, 42, 43, 44, 45, 46, 47 48, 49, 50, 51, 52, 53, 54, 55

PROSTHESIS

ACL SCREW



►DIA Ø ►REF. NO. TT

7.0 MM 26-7020T to 26-7040T

►LENGTH IN MM 20, 25, 30, 35, 40

8.0 MM 26-8020T to 26-8040T

►LENGTH IN MM 20, 25, 30, 35, 40

ACL SCREW



►DIA Ø ►REF. NO. TT

9.0 MM 26-9020T to 26-9040T

►LENGTH IN MM 20, 25, 30, 35, 40

10 MM 26-1020T to 26-1040T

►LENGTH IN MM 20, 25, 30, 35, 40

ACL ENDO BUTTON



►REF. NO. TT

26-0001

ACL DISC WASHER



►REF. NO. TT

26-0002

LOCKING PLATES



LOCKING PLATES

3.5 MM LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION

RIGHT SIDE

►REF. NO. SS	13-3506 to 13-3508
►REF. NO. TT	13-3506T to 13-3508T

LEFT SIDE

►REF. NO. SS	13-3506 to 13-3508
►REF. NO. TT	13-3506T to 13-3508T
►HOLES	6, 7, 8



3.5 MM LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION

RIGHT SIDE

►REF. NO. SS	13-3503 to 13-3508
►REF. NO. TT	13-3503T to 13-3508T

LEFT SIDE

►REF. NO. SS	13-3503 to 13-3508
►REF. NO. TT	13-3503T to 13-3508T
►HOLES	3, 4, 5, 6, 7, 8



3.5 MM LCP CLAVICLE HOOK PLATE

RIGHT SIDE

►REF. NO. SS	13-3504 to 13-3507
►DEPTH	12 MM, 15 MM, 18 MM

LEFT SIDE

►REF. NO. SS	13-3504 to 13-3507
►DEPTH	12 MM, 15 MM, 18 MM
►HOLES	4, 5, 6, 7



3.5 MM LCP SMALL PLATE

►REF. NO. SS	13-3505 to 13-3514
►REF. NO. TT	13-3505T to 13-3514T
►HOLES	5, 6, 7, 8, 9, 10, 11, 12, 14



LOCKING PLATE

3.5 MM LCP PERIARTICULAR PROXIMAL HUMERUS PLATE



RIGHT SIDE

- ▶REF. NO. SS 13-3502 to 13-3514
- ▶REF. NO. TT 13-3502T to 13-3514T

LEFT SIDE

- ▶REF. NO. SS 13-3502 to 13-3514
- ▶REF. NO. TT 13-3502T to 13-3514T
- ▶HOLES 2, 4, 6, 8, 10, 12, 14

3.5 MM LONG PHILOS PLATE



- ▶REF. NO. SS 13-3503 to 13-3513
- ▶REF. NO. TT 13-3503T to 13-3513T
- ▶HOLES 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13

3.5 MM SHORT PHILOS PLATE



- ▶REF. NO. SS 13-3503 to 13-3505
- ▶REF. NO. TT 13-3503T to 13-3505T
- ▶HOLES 3, 5

3.5 MM LCP RECONSTRUCTION PLATE



- ▶REF. NO. SS 13-3504 to 13-3522
- ▶REF. NO. TT 13-3504T to 13-3522T
- ▶HOLES 4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18, 20, 22

LOCKING PLATES

3.5 / 5 MM LCP METAPHYSEAL PLATE

►REF. NO. SS	13-3508 to 13-3518
►REF. NO. TT	13-3508T to 13-3518T
►HOLES	3, 4, 5, 6, 7, 8, 9, 11, 13

3.5 MM OLECRANON PLATE

RIGHT SIDE

►REF. NO. SS	13-3504 to 13-3512
►REF. NO. TT	13-3504T to 13-3512T

LEFT SIDE

►REF. NO. SS	13-3504 to 13-3512
►REF. NO. TT	13-3504T to 13-3512T
►HOLES	4, 6, 8, 10, 12

3.5 MM LCP MEDIAL DISTAL HUMERUS PLATE

RIGHT SIDE

►REF. NO. SS	13-3503 to 13-3514
►REF. NO. TT	13-3503T to 13-3514T

LEFT SIDE

►REF. NO. SS	13-3503 to 13-3514
►REF. NO. TT	13-3503T to 13-3514T
►HOLES	3, 5, 7, 9, 14

3.5 MM LCP POSTEROLATERAL DISTAL HUMERUS PLATE WITH LATERAL SUPPORT

RIGHT SIDE

►REF. NO. SS	13-3503 to 13-3514
►REF. NO. TT	13-3503T to 13-3514T

LEFT SIDE

►REF. NO. SS	13-3503 to 13-3514
►REF. NO. TT	13-3503T to 13-3514T
►HOLES	3, 5, 7, 9, 14

LOCKING PLATES

3.5 MM LCP POSTEROLATERAL DISTAL HUMERUS PLATE WITHOUT LATERAL SUPPORT

RIGHT SIDE

►REF. NO. SS	13-3503 to 13-3514
►REF. NO. TT	13-3503T to 13-3514T

LEFT SIDE

►REF. NO. SS	13-3503 to 13-3514
►REF. NO. TT	13-3503T to 13-3514T
►HOLES	3, 5, 7, 9, 14



3.5 MM LCP EXTRA - ARTICULAR DISTAL HUMERUS PLATE

RIGHT SIDE

►REF. NO. SS	13-3504 to 13-3514
►REF. NO. TT	13-3504T to 13-3514T

LEFT SIDE

►REF. NO. SS	13-3504 to 13-3514
►REF. NO. TT	13-3504T to 13-3514T
►HOLES	4, 6, 8, 10, 12, 14



2.4 MM / 2.7 MM VA - TWO COLUMN VOLAR DISTAL RADIUS PLATE

►REF. NO. SS	13-2702 to 13-2705
►REF. NO. TT	13-2702T to 13-2705T
►HOLES	2, 3, 4, 5



3.5 MM LCP T-PLATE RIGHT ANGLE

►REF. NO. SS	13-3503 to 13-3508
►REF. NO. TT	13-3503T to 13-3508T
►HOLES	3, 4, 5, 6, 7, 8



LOCKING PLATES

3.5 MM LCP T-PLATE OBLIQUE ANGLE

RIGHT SIDE

- ▶REF. NO. SS 13-3503 to 13-3508
- ▶REF. NO. TT 13-3503T to 13-3508T

LEFT SIDE

- ▶REF. NO. SS 13-3503 to 13-3508
- ▶REF. NO. TT 13-3503T to 13-3508T
- ▶HOLES 3, 4, 5, 6, 7, 8

2.4 MM LCP DISTAL RADIUS PLATE STRAIGHT - DORSAL

- ▶REF. NO. SS 13-2405 to 13-2406
- ▶REF. NO. TT 13-2405T to 13-2406T
- ▶HOLES 5, 6

2.4 MM LCP L DISTAL RADIUS PLATE - DORSAL

RIGHT SIDE

- ▶REF. NO. SS 13-2403 to 13-2404
- ▶REF. NO. TT 13-2403T to 13-2404T

LEFT SIDE

- ▶REF. NO. SS 13-2403 to 13-2404
- ▶REF. NO. TT 13-2403T to 13-2404T
- ▶HOLES 2 HEAD HOLES, 3, 4 SHAFT HOLES

2.4 MM LCP L DISTAL RADIUS PLATE - DORSAL

RIGHT SIDE

- ▶REF. NO. SS 13-2403 to 13-2404
- ▶REF. NO. TT 13-2403T to 13-2404T

LEFT SIDE

- ▶REF. NO. SS 13-2403 to 13-2404
- ▶REF. NO. TT 13-2403T to 13-2404T
- ▶HOLES 3 HEAD HOLES, 3, 4 SHAFT HOLES

LOCKING PLATES

2.4 MM LCP L DISTAL RADIUS PLATE OBLIQUE - DORSAL



RIGHT SIDE

- ▶REF. NO. SS 13-2403 to 13-2404
- ▶REF. NO. TT 13-2403T to 13-2404T

LEFT SIDE

- ▶REF. NO. SS 13-2403 to 13-2404
- ▶REF. NO. TT 13-2403T to 13-2404T
- ▶HOLES 3 HEAD HOLES, 3, 4 SHAFT HOLES

2.4 MM LCP T DISTAL RADIUS PLATE - DORSAL



- ▶REF. NO. SS 13-2403 to 13-2404
- ▶REF. NO. TT 13-2403T to 13-2404T
- ▶HOLES 3 HEAD HOLES, 3, 4 SHAFT HOLES

2.4 MM LCP DISTAL RADIUS PLATE - VOLAR



RIGHT SIDE

- ▶REF. NO. SS 13-2403 to 13-2405
- ▶REF. NO. TT 13-2403T to 13-2405T

LEFT SIDE

- ▶REF. NO. SS 13-2403 to 13-2405
- ▶REF. NO. TT 13-2403T to 13-2405T
- ▶HOLES 5 HEAD HOLES - 3 SHAFT HOLES
5 HEAD HOLES - 5 SHAFT HOLES

2.4 MM LCP BUTTRESS PLATE - VOLAR



RIGHT SIDE

- ▶REF. NO. SS 13-2403
- ▶REF. NO. TT 13-2403T

LEFT SIDE

- ▶REF. NO. SS 13-2403
- ▶REF. NO. TT 13-2403T
- ▶HOLES 5 HEAD HOLES - 3 SHAFT HOLES

LOCKING PLATES

2.4 MM LCP DISTAL RADIUS - EXTRAARTICULAR PLATE - VOLAR



RIGHT SIDE

- ▶REF. NO. SS 13-2403 to 13-2405
- ▶REF. NO. TT 13-2403T to 13-2405T

LEFT SIDE

- ▶REF. NO. SS 13-2403 to 13-2405
- ▶REF. NO. TT 13-2403T to 13-2405T
- ▶HOLES 5 HEAD HOLES, 3, 5 SHAFT HOLES

2.4 MM LCP DISTAL RADIUS - EXTRAARTICULAR PLATE- VOLAR



RIGHT SIDE

- ▶REF. NO. SS 13-2403 to 13-2405
- ▶REF. NO. TT 13-2403T to 13-2405T

LEFT SIDE

- ▶REF. NO. SS 13-2403 to 13-2405
- ▶REF. NO. TT 13-2403T to 13-2405T
- ▶HOLES 4 HEAD HOLES, 3, 5 SHAFT HOLES

2.4 MM LCP DISTAL RADIUS - EXTRA LONG PLATE - VOLAR



- ▶REF. NO. SS 13-2408 to 13-2412
- ▶REF. NO. TT 13-2408T to 13-2412T
- ▶HOLES 4 HEAD HOLES - 8, 10 & 12 SHAFT HOLES

5.0 MM LCP BROAD PLATE



- ▶REF. NO. SS 13-5005 to 13-5018
- ▶REF. NO. TT 13-5005T to 13-5018T
- ▶HOLES 5, 6, 7, 8, 9, 10, 11, 12, 14, 16, 18

LOCKING PLATES

5.0 MM LCP DISTAL FEMUR PLATE



RIGHT SIDE

- ▶REF. NO. SS 13-5005 to 13-5015
- ▶REF. NO. TT 13-5005T to 13-5015T

LEFT SIDE

- ▶REF. NO. SS 13-5005 to 13-5015
- ▶REF. NO. TT 13-5005T to 13-5015T
- ▶HOLES 5, 7, 9, 11, 13, 15

5.0 MM LCP PROXIMAL FEMUR PLATE (TROCHANTRIC PLATE)



- ▶REF. NO. SS 13-5005 to 13-5013
- ▶REF. NO. TT N/A
- ▶HOLES 5, 7, 9, 11, 13

5.0 MM LCP NARROW PLATE



- ▶REF. NO. SS 13-5005 to 13-5016
- ▶REF. NO. TT 13-5005T to 13-5016T
- ▶HOLES 5, 6, 7, 8, 9, 10, 11, 12, 14, 16

5.0 MM LCP T-BUTTRESS PLATE



- ▶REF. NO. SS 13-5004 to 13-5012
- ▶REF. NO. TT 13-5004T to 13-5012T
- ▶HOLES 4, 5, 6, 7, 8, 9, 10, 12

LOCKING PLATES

5.0 MM LCP L- BUTTRESS PLATE

RIGHT SIDE

- ▶REF. NO. SS 13-5004 to 13-5012
- ▶REF. NO. TT 13-5004T to 13-5012T

LEFT SIDE

- ▶REF. NO. SS 13-5004 to 13-5012
- ▶REF. NO. TT 13-5004T to 13-5012T
- ▶HOLES 4, 5, 6, 7, 8, 9, 10, 12

3.5 MM LCP LOW BEND MEDIAL DISTAL TIBIA PLATE

RIGHT SIDE

- ▶REF. NO. SS 13-3504 to 13-3514
- ▶REF. NO. TT 13-3504T to 13-3514T

LEFT SIDE

- ▶REF. NO. SS 13-3504 to 13-3514
- ▶REF. NO. TT 13-3504T to 13-3514T
- ▶HOLES 4, 6, 8, 10, 12, 14

3.5 MM LCP MEDIAL DISTAL TIBIA PLATE - LEFT SIDE (WITH TAB)

RIGHT SIDE

- ▶REF. NO. SS 13-3504 to 13-3514
- ▶REF. NO. TT 13-3504T to 13-3514T

LEFT SIDE

- ▶REF. NO. SS 13-3504 to 13-3514
- ▶REF. NO. TT 13-3504T to 13-3514T
- ▶HOLES 4, 6, 8, 10, 12, 14

5.0 MM LCP PROXIMAL LATERAL TIBIA PLATE - (HOCKEY PLATE)

RIGHT SIDE

- ▶REF. NO. SS 13-5005 to 13-5015
- ▶REF. NO. TT 13-5005T to 13-5015T

LEFT SIDE

- ▶REF. NO. SS 13-5005 to 13-5015
- ▶REF. NO. TT 13-5005T to 13-5015T
- ▶HOLES 5, 7, 9, 11, 13, 15

LOCKING PLATES

5.0 MM LCP RECONSTRUCTION PLATE



►REF. NO. SS	13-5004 to 13-5016
►REF. NO. TT	13-5004T to 13-5016T
►HOLES	4, 5, 6, 7, 8, 9, 10, 12, 14, 16

5.0 MM LCP T-PLATE



►REF. NO. SS	13-5003 to 13-5012
►REF. NO. TT	13-5003T to 13-5012T
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

3.5 MM LCP CLOVERLEAF PLATE



►REF. NO. SS	13-3503 to 13-3512
►REF. NO. TT	13-3503T to 13-3512T
►HOLES	3, 4, 5, 6, 7, 8, 9, 10, 12

3.5 MM LCP LATERAL DISTAL FIBULA PLATE



►REF. NO. SS	13-3503 to 13-3511
►REF. NO. TT	13-3503T to 13-3511T
►HOLES	3, 4, 5, 6, 7, 9, 10, 11

LOCKING PLATES

3.5 MM LCP PROXIMAL TIBIA PLATE (RAFT)



RIGHT SIDE

- ▶REF. NO. SS 13-3504 to 13-3516
- ▶REF. NO. TT 13-3504T to 13-3516T

LEFT SIDE

- ▶REF. NO. SS 13-3504 to 13-3516
- ▶REF. NO. TT 13-3504T to 13-3516T
- ▶HOLES 4, 6, 8, 10, 12, 14, 16

3.5 MM ANTEROLATERAL DISTAL TIBIA PLATE



RIGHT SIDE

- ▶REF. NO. SS 13-3505 to 13-3517
- ▶REF. NO. TT 13-3505T to 13-3517T

LEFT SIDE

- ▶REF. NO. SS 13-3505 to 13-3517
- ▶REF. NO. TT 13-3505T to 13-3517T
- ▶HOLES 5, 7, 9, 11, 13, 15, 17

3.5 MM LCP MEDIAL PROXIMAL TIBIA PLATE



RIGHT SIDE

- ▶REF. NO. SS 13-3504 to 13-3516
- ▶REF. NO. TT 13-3504T to 13-3516T

LEFT SIDE

- ▶REF. NO. SS 13-3504 to 13-3516
- ▶REF. NO. TT 13-3504T to 13-3516T
- ▶HOLES 4, 6, 8, 10, 12, 14, 16

3.5 MM LCP CALCANEAL PLATE



RIGHT SIDE

- ▶REF. NO. SS 13-3564 to 13-3581
- ▶REF. NO. TT 13-3564T to 13-3581T

LEFT SIDE

- ▶REF. NO. SS 13-3564 to 13-3581
- ▶REF. NO. TT 13-3564T to 13-3581T
- ▶HOLES 64, 69, 76, 81

MAXILLOFACIAL IMPLANTS



1.5 MM STRAIGHT PLATE WITH STEM/ BAR/ GAP - 0.6 MM THICKNESS



►REF. NO. SS	N/A
►REF. NO. TT	24-1502T to 24-1506T
►HOLES	2, 4, 6

2.0 MM STRAIGHT PLATE RECON WITH STEM/ BAR/ GAP - 1 MM THICKNESS



►REF. NO. SS	N/A
►REF. NO. TT	24-2002 to 24-2006T
►HOLES	2, 4, 6

1.5 MM STRAIGHT PLATE CONTINUES RECON - 0.6 MM THICKNESS



►REF. NO. SS	N/A
►REF. NO. TT	24-1502T to 24-1518T
►HOLES	2, 4, 6, 18

2.0 MM STRAIGHT PLATE CONTINUES RECON - 1 MM THICKNESS



►REF. NO. SS	N/A
►REF. NO. TT	24-2002T to 24-2018T
►HOLES	2, 4, 6, 18

MAXILLOFACIAL

1.5 MM T - PLATE RECON - 0.6 MM THICKNESS



►REF. NO. SS	N/A
►REF. NO. TT	24-1502T to 24-1503T
►HOLES	3 HEAD HOLES + 2 SHAFT HOLES
	3 HEAD HOLES + 3 SHAFT HOLES

2.0 MM T - PLATE RECON - 1.0 MM THICKNESS



►REF. NO. SS	N/A
►REF. NO. TT	24-2002T to 24-2003T
►HOLES	3 HEAD HOLES + 2 SHAFT HOLES
	3 HEAD HOLES + 3 SHAFT HOLES

1.5 MM T-OBLIQUE PLATE RECON WITH GAP - 0.6 MM THICKNESS

►REF. NO. SS	N/A
►REF. NO. TT	24-1502T
►HOLES	3 HEAD HOLES + 2 SHAFT HOLES

2.0 MM T-OBLIQUE PLATE RECON WITH GAP - 1.0 MM THICKNESS

►REF. NO. SS	N/A
►REF. NO. TT	24-2002T
►HOLES	3 HEAD HOLES + 2 SHAFT HOLES

MAXILLOFACIAL

1.5 MM DOUBLE T-PLATE RECON - 0.6 MM THICKNESS

▶REF. NO. SS	N/A
▶REF. NO. TT	24-1508T
▶HOLES	8

2.0 MM DOUBLE T-PLATE RECON - 1.0 MM THICKNESS

▶REF. NO. SS	N/A
▶REF. NO. TT	24-2008T
▶HOLES	8

1.5 MM DOUBLE SQUARE PLATE RECON - 0.6 MM THICKNESS

▶REF. NO. SS	N/A
▶REF. NO. TT	24-1506T
▶HOLES	6

2.0 MM DOUBLE SQUARE PLATE RECON - 1.0 MM THICKNESS

▶REF. NO. SS	N/A
▶REF. NO. TT	24-2006T
▶HOLES	6




MAXILLOFACIAL

1.5 MM DOUBLE Y- PLATE RECON

WITH STEM/ BAR/ GAP - 0.6 MM THICKNESS


►REF. NO. SS	N/A
►REF. NO. TT	24-1506T
►HOLES	6



2.0 MM DOUBLE Y- PLATE RECON

WITH STEM/ BAR/ GAP - 1.0 MM THICKNESS

►REF. NO. SS	N/A
►REF. NO. TT	24-2006T
►HOLES	6



1.5 MM L-PLATE RECON - 0.6 MM THICKNESS

RIGHT SIDE

►REF. NO. TT	24-1502T to 24-1503T
►HOLES	3 HEAD HOLES + 2 SHAFT HOLES 3 HEAD HOLES + 3 SHAFT HOLES

LEFT SIDE

►REF. NO. TT	24-1502T to 24-1503T
►HOLES	3 HEAD HOLES + 2 SHAFT HOLES 3 HEAD HOLES + 3 SHAFT HOLES

2.0 MM L-PLATE RECON - 1.0 MM THICKNESS - RIGHT SIDE

RIGHT SIDE

►REF. NO. TT	24-2002T to 24-2003T
►HOLES	3 HEAD HOLES + 2 SHAFT HOLES 3 HEAD HOLES + 3 SHAFT HOLES

LEFT SIDE

►REF. NO. TT	24-2002T to 24-2003T
►HOLES	3 HEAD HOLES + 2 SHAFT HOLES 3 HEAD HOLES + 3 SHAFT HOLES



MAXILLOFACIAL

1.5 MM L-OBLIQUE PLATE RECON - 0.6 THICKNESS - RIGHT

RIGHT SIDE

- ▶REF. NO. TT 24-1503T
- ▶HOLES 2 HEAD HOLES + 3 SHAFT HOLES

LEFT SIDE

- ▶REF. NO. TT 24-1503T
- ▶HOLES 2 HEAD HOLES + 3 SHAFT HOLES



2.0 MM L-OBLIQUE PLATE RECON - 1.0 THICKNESS - RIGHT SIDE - 120°

RIGHT SIDE

- ▶REF. NO. TT 24-2003T
- ▶HOLES 2 HEAD HOLES + 3 SHAFT HOLES

LEFT SIDE

- ▶REF. NO. TT 24-2003T
- ▶HOLES 2 HEAD HOLES + 3 SHAFT HOLES



1.5 MM SQUARE PLATE RECON - 0.6 MM THICKNESS

- ▶REF. NO. SS N/A
- ▶REF. NO. TT 24-1504T
- ▶HOLES 4




2.0 MM SQUARE PLATE RECON - 1.0 MM THICKNESS

- ▶REF. NO. SS N/A
- ▶REF. NO. TT 24-2004T
- ▶HOLES 4

MAXILLOFACIAL

1.5 MM X - PLATE RECON - 0.6 MM THICKNESS



►REF. NO. SS	N/A
►REF. NO. TT	2-1502T to 24-1504T
►HOLES	2 HEAD HOLES + 2 SHAFT HOLES
	4 HEAD HOLES + 4 SHAFT HOLES

2.0 MM X - PLATE RECON - 1.0 MM THICKNESS



►REF. NO. SS	N/A
►REF. NO. TT	24-2002T to 24-2004T
►HOLES	2 HEAD HOLES + 2 SHAFT HOLES
	4 HEAD HOLES + 4 SHAFT HOLES

1.5 MM Y- PLATE RECON



►REF. NO. SS	N/A
►REF. NO. TT	24-1500T to 24-1500T
►HOLES	SMALL
	BIG

2.0 MM Y- PLATE RECON




►REF. NO. SS	N/A
►REF. NO. TT	24-2000T to 24-2000T
►HOLES	SMALL
	BIG



MAXILLOFACIAL

1.5 MM ORBITAL PLATE RECON

►REF. NO. SS	N/A
►REF. NO. TT	24-1504T to 24-1507T
►HOLES	4H (2+2)
	6H (3+3)
	7H (4+3)



2.0 MM ORBITAL PLATE RECON CONTINUES

►REF. NO. SS	N/A
►REF. NO. TT	24-2008T
►HOLES	8



1.5 MM MAXILLOFACIAL SCREWS (SELF - TAPPING) - HEAD CROSS SLOT

►REF. NO. SS	N/A
►REF. NO. TT	24-1505T to 24-1516T
►HOLES	5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16



2.0 MM MAXILLOFACIAL SCREWS (SELF - TAPPING) - HEAD CROSS SLOT

►REF. NO. SS	N/A
►REF. NO. TT	24-2005T to 24-2034T
►HOLES	5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17
	18, 19, 20, 21, 22, 23, 24, 25, 26, 27, 28, 29
	30, 31, 32, 33, 34

PELVIC IMPLANTS



PELVIC IMPLANTS

3.5 MM PELVIC RECON PLATE LOW PROFILE - STRAIGHT



►REF. NO. SS	22-3502 to 22-3520
►REF. NO. TT	N/A
►HOLES	2 MM to 10 MM WITH 1 MM DIFF. 10 MM to 20 MM WITH 2 MM DIFF.

3.5 MM PELVIC RECON PLATE 88° - CURVED



►REF. NO. SS	22-3506 to 22-3516
►REF. NO. TT	N/A
►HOLES	6, 8, 10, 12, 14, 16

3.5 MM PELVIC RECON PLATE 108° - CURVED



►REF. NO. SS	22-3506 to 22-3516
►REF. NO. TT	N/A
►HOLES	6, 8, 10, 12, 14, 16

3.5 MM LCP PELVIC RECON PLATE - STRAIGHT



►REF. NO. SS	22-3502 to 22-3518
►REF. NO. TT	N/A
►HOLES	2 MM to 10 MM WITH 1 MM DIFF. 10 MM to 18 MM WITH 2 MM DIFF.

PELVIC IMPLANTS

3.5 MM SPRING PLATE



►REF. NO. SS	22-3501 to 22-3503
►REF. NO. TT	N/A
►HOLES	1 HOLE - LENGTH - 19.5 MM
	2 HOLES - LENGTH - 31.5 MM
	3 HOLES - LENGTH - 43.5 MM

3.5 MM PELVIC J PLATE



RIGHT SIDE

►REF. NO. SS	22-3510 to 22-3516
►REF. NO. TT	N/A
►HOLES	10, 12, 14, 16

LEFT SIDE

►REF. NO. SS	22-3510 to 22-3516
►REF. NO. TT	N/A
►HOLES	10, 12, 14, 16

PINS & WIRES IMPLANTS



PINS & WIRES

3.5 MM MOORE'S PIN WITH 2 NUTS



- REF. NO. SS 30-3045 to 30-30150
- LENGTH IN MM 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95
100, 105, 110, 115, 120, 125, 130, 135
140, 145, 150

KNOWLES PIN



- DIA Ø ►REF. NO. SS
3.0 MM 30-30100 to 30-30145
- LENGTH IN MM 100, 105, 110, 115, 120, 125, 130, 135
140, 145
- 4.0 MM 30-4095 to 30-40180
- LENGTH IN MM 95, 100, 105, 110, 115, 120, 125, 130, 135
140, 145, 150, 155, 160, 165, 170, 175, 180

2.5 MM SCHANZ PIN



- DIA Ø ►REF. NO. SS
2.5 MM 30-14100
- 14 THREAD 100 IN MM

3.5 MM SCHANZ PIN



- DIA Ø ►REF. NO. SS
3.5 MM 30-16125 to 30-40225
- 16 THREAD 125, 150, 175, 200 IN MM.
- 18 THREAD 175 IN MM
- 32 THREAD 150, 175, 200, 225 IN MM.
- 40 THREAD 225 IN MM.

PINS & WIRES

3.5 MM X 4.0 MM SHAFT SCHANZ PIN

►DIA Ø ►REF. NO. SS

3.5 MM 30-18175 to 30-32150

►18 THREAD 175 IN MM.

►25 THREAD 125 IN MM.

►32 THREAD 150 IN MM.

3.5 MM X 4.5 MM SHAFT SCHANZ PIN

►REF. NO. SS 30-32125 to 30-32150

►32 THREAD 125, 150 IN MM.

4.0 MM SCHANZ PIN

►REF. NO. SS 30-40175

►40 THREAD 175 IN MM.

4.5 MM SCHANZ PIN

►DIA Ø ►REF. NO. SS

4.5 MM 30-18150 to 30-40225

►18 THREAD 150, 175, 200, 225 IN MM.

►25 THREAD 125, 150, 175, 200, 225, 250 IN MM

►32 THREAD 150, 175, 200, 225 IN MM.

►40 THREAD 175, 200, 225 IN MM.

PINS & WIRES

5.0 MM SCHANZ PIN

►DIA Ø ►REF. NO. SS

5.0 MM 30-16175 to 30-50225

►16 THREAD 175 IN MM.

►18 THREAD 150, 175, 200 IN MM.

►25 THREAD 150, 175, 200 IN MM.

►32 THREAD 150, 175, 200, 225 IN MM.

►50 THREAD 100, 125, 150, 175, 200, 225 IN MM.

6.0 MM SCHANZ PIN

►DIA Ø ►REF. NO. SS

6.0 MM 30-32150 to 30-50225

►32 THREAD 150, 175, 200, 225 IN MM.

►50 THREAD 125, 150, 175, 200, 225 IN MM.

4.0 MM TAPERED SCHANZ PIN

►DIA Ø ►REF. NO. SS

4.0 MM 30-4050 to 30-40100

►18 THREAD 50 IN MM.

►20 THREAD 60, 70, 80, 100 IN MM.

4.5 MM TAPERED SCHANZ PIN

►DIA Ø ►REF. NO. SS

4.5 MM 30-20100 to 30-40170

►20 THREAD 100 IN MM.

►30 THREAD 100, 120 IN MM.

►40 THREAD 120, 150, 170 IN MM.

PINS & WIRES

6.0 MM TAPERED SCHANZ PIN

►DIA Ø ►REF. NO. SS

6.0 MM 30-30130 to 30-60200

►30 THREAD 130, 150 IN MM.

►40 THREAD 130, 150, 170, 200 IN MM.

►50 THREAD 130, 150, 170, 200, 225, 250 IN MM.

►60 THREAD 170, 200 IN MM.

SCHANZ PIN SELF DRILLING AND SELF TAPPING

►REF. NO. SS

30-40150 to 30-75200

►40 THREAD 150 IN MM.

►75 THREAD 195, 200 IN MM.

SCHANZ PIN (STERILISED)

►REF. NO. SS 30-16150 to 30-25200

►16 THREAD 150 IN MM.

►18 THREAD 200 IN MM.

►25 THREAD 200 IN MM.

6.5 MM CANCELLOUS SCHANZ PIN

►DIA Ø ►REF. NO. SS

6.5 MM 30-16225 to 30-32225

►16 THREAD 225 IN MM.

►32 THREAD 225 IN MM.

PINS & WIRES

CANCELLOUS SCHANZ PIN WITH THREADED END WITH 2 NUTS 6.5 MM

►REF. NO. SS

30-0001

►LENGTH IN MM 150



STEINMANN PIN

►DIA Ø ►REF. NO. SS

3.0 MM 30-30150 to 30-30225

►LENGTH IN MM 150, 175, 200, 225

3.5 MM 30-35125 to 30-35250

►LENGTH IN MM 125, 150, 175, 200, 225, 250

4.0 MM 30-40150 to 30-40250

►LENGTH IN MM 150, 175, 200, 225, 250



STEINMANN PIN

►DIA Ø ►REF. NO. SS

4.5 MM 30-45150 to 30-45250

►LENGTH IN MM 150, 175, 200, 225, 250

5.0 MM 30-50150 to 30-50225

►LENGTH IN MM 150, 175, 200, 225



4.5 MM STEINMANN PIN (STERILISED)

►DIA Ø ►REF. NO. SS

4.5 MM 30-45225

►LENGTH IN MM 225



PINS & WIRES

DANHAM PIN

►DIA Ø ►REF. NO. SS

4.5 MM 30-45225

►LENGTH IN MM 225

5.0 MM 30-50225

►LENGTH IN MM 225

DANHAM PIN (STERILIZED)

►DIA Ø ►REF. NO. SS

4.5 MM 30-45225

►LENGTH IN MM 225

K. WIRE

►DIA Ø ►REF. NO. SS

0.6 MM 30-06100 to 30-06250

►LENGTH IN MM 100, 150, 200, 250

0.8 MM 30-08075 to 30-08150

►LENGTH IN MM 75, 100, 150

1.0 MM 30-10100 to 30-10300

►LENGTH IN MM 100, 150, 200, 250, 300

K. WIRE

►DIA Ø ►REF. NO. SS

1.2 MM 30-12100 to 30-12300

►LENGTH IN MM 100, 150, 200, 250, 300

1.4 MM 30-14100 to 30-14310

►LENGTH IN MM 100, 150, 200, 250, 300, 310

1.5 MM 30-15100 to 30-15400

►LENGTH IN MM 100, 150, 200, 250, 300, 350, 400

PINS & WIRES

K. WIRE

- ▶DIA Ø ▶REF. NO. SS
- 1.6 MM 30-16100 to 30-16310
- ▶LENGTH IN MM 100, 150, 200, 250, 300, 310
- 1.7 MM 30-17150
- ▶LENGTH IN MM 150
- 1.8 MM 30-18100 to 30-18310
- ▶LENGTH IN MM 100, 150, 175, 200, 250, 300, 310

K. WIRE

- ▶DIA Ø ▶REF. NO. SS
- 2.0 MM 30-20100 to 30-20400
- ▶LENGTH IN MM 100, 150, 175, 200, 225, 250, 300, 310, 400
- 2.5 MM 30-25100 to 30-25450
- ▶LENGTH IN MM 100, 150, 200, 250, 300, 400, 450
- 3.0 MM 30-30100 to 30-30450
- ▶LENGTH IN MM 100, 150, 200, 250, 300, 310, 450

K. WIRE

- ▶DIA Ø ▶REF. NO. SS
- 3.5 MM 30-35100 to 30-35300
- ▶LENGTH IN MM 100, 150, 200, 250, 300
- 4.0 MM 30-40150 to 30-40350
- ▶LENGTH IN MM 150, 200, 250, 300, 350
- 5.0 MM 30-50300
- ▶LENGTH IN MM 300

K. WIRE (STERILIZED)

- ▶DIA Ø ▶REF. NO. SS
- 1.0 MM 30-10150
- ▶LENGTH IN MM 150
- 1.2 MM 30-12150
- ▶LENGTH IN MM 150
- 1.5 MM 30-15150
- ▶LENGTH IN MM 150

PINS & WIRES

K. WIRE (STERILIZED)

►DIA Ø ►REF. NO. SS

1.8 MM 30-18150

►LENGTH IN MM 150

2.0 MM 30-20150

►LENGTH IN MM 150

2.5 MM 30-25150

►LENGTH IN MM 150

K. WIRE (STERILIZED)

►DIA Ø ►REF. NO. SS

3.0 MM 30-30150

►LENGTH IN MM 150

DOUBLE ENDED- THREADED K. WIRE

►DIA Ø ►REF. NO. SS

1.2 MM 30-12300

►LENGTH IN MM 300

1.4 MM 30-14250

►LENGTH IN MM 250

1.5 MM 30-15150 to 30-15300

►LENGTH IN MM 150, 180, 200, 250, 300

DOUBLE ENDED- THREADED K. WIRE

►DIA Ø ►REF. NO. SS

2.0 MM 30-20150 to 30-20250

►LENGTH IN MM 150, 200, 250

2.5 MM 30-25100 to 30-25400

►LENGTH IN MM 100, 150, 200, 250, 300, 400

3.0 MM 30-30150 to 30-30350

►LENGTH IN MM 150, 200, 250, 300, 350

PINS & WIRES

3.5 MM DOUBLE ENDED- THREADED K. WIRE



►DIA Ø ►REF. NO. SS

3.5 MM 30-35300

►LENGTH IN MM 300

S.S. WIRE (SUTURE)



►REF. NO. SS

30-0016 to 30-0030

►LENGTH IN MM 16, 18, 20, 22, 24, 26, 28, 30

COVENTRY STAPLES



►REF. NO. SS

30-0500 to 30-1500

►SIZE: SMALL 5 MM STEP

MEDIUM 10 MM STEP

LARGE 15 MM STEP

PINS & WIRES

ILLIZAROV WIRE 1.5 (BAYNOT & TROCAR POINT,

BAYNOT POINT



TROCAR POINT



►DIA Ø ►REF. NO. SS

1.5 MM 30-1512, 30-1516

►LENGTH IN MM 12 & 16

1.8 MM 30-1816

►LENGTH IN MM 16

OLIVE WIRE 1.5 (BAYNOT & TROCAR POINT)

BAYNOT POINT



TROCAR POINT



►DIA Ø ►REF. NO. SS

1.5 MM 30-1516

►BAYNOT POINT 16

1.8 MM 30-1816

►BAYNOT POINT 16

SPINAL IMPLANTS



SPINAL IMPLANTS

MONOAXIAL SCREW



►DIA Ø ►REF. NO. TT

4.5 MM 20-4525T to 20-4550T

►LENGTH IN MM 25, 30, 35, 40, 45, 50

5.5 MM 20-5525T to 20-5550T

►LENGTH IN MM 25, 30, 35, 40, 45, 50

6.5 MM 20-6525T to 20-6550T

►LENGTH IN MM 25, 30, 35, 40, 45, 50

POLYAXIAL SCREW



►DIA Ø ►REF. NO. TT

4.5 MM 20-4525T to 20-4550T

►LENGTH IN MM 25, 30, 35, 40, 45, 50

5.5 MM 20-5530T to 20-5550T

►LENGTH IN MM 30, 35, 40, 45, 50

6.5 MM 20-6530T to 20-6550T

►LENGTH IN MM 30, 35, 40, 45, 50

REDUCTION POLYAXIAL SCREW



►DIA Ø ►REF. NO. TT

5.5 MM 20-5530T to 20-5550T

►LENGTH IN MM 30, 35, 40, 45, 50

6.5 MM 20-6530T to 20-6550T

►LENGTH IN MM 30, 35, 40, 45, 50

LATERAL MASS SCREW (POSTERIOR CERVICAL POLY)



►DIA Ø ►REF. NO. TT

3.5 MM 20-3510T to 20-3530T

4.5 MM 20-4510T to 20-4530T

►LENGTH IN MM 10 MM to 30 MM WITH 2 MM DIFF.

SPINAL IMPLANTS

CERVICAL PLATE



►REF. NO. SS

20-0022T to 20-0070T

►LENGTH IN MM 22.5, 25, 27.5, 30, 32.2, 35, 37.5, 40

42.5, 45, 47.5, 50, 55, 60, 65, 70

CERVICAL SCREW



►DIA Ø

3.5 MM

►REF. NO. TT

20-3510T to 20-3520T

►LENGTH IN MM 10, 12, 14, 16, 18, 20

4.0 MM

2025-TT-4010 to 2025-TT-4020

►LENGTH IN MM 10, 12, 14, 16, 18, 20

CERVICAL PLATE WITH INBUILT CAGE



►THICKNESS

5.0 MM

►REF. NO. TT

20-5000T

6.0 MM

20-6000T

7.0 MM

20-7000T

8.0 MM

20-8000T

CERVICAL PLATE WITH INBUILT CAGE SCREW



►DIA Ø

3.5 MM

►REF. NO. TT

20-3510T to 20-3520T

►LENGTH IN MM

10, 12, 14, 16, 18, 20

3.8 MM

20-3810T to 20-3820T

►LENGTH IN MM

10, 12, 14, 16, 18, 20

SPINAL IMPLANTS

SPINAL ROD FOR LATERAL MASS.

►DIA Ø ►REF. NO. TT

3.5 MM 20-3540T to 20-35250T

►LENGTH IN MM 40 MM to 200 MM WITH 10 MM DIFF.
200, 250



SPINAL ROD FOR DORSOLUMBER

►DIA Ø ►REF. NO. TT

5.5 MM 20-5540T to 20-55500T

►LENGTH IN MM 40 MM to 300 MM WITH 10 MM DIFF.
300 MM to 500 MM WITH 50 MM DIFF.



ECO CAGE

►DIA Ø ►REF. NO. TT

10 MM 20-1010T to 20-1050T

12 MM 20-1210T to 20-1250T

14 MM 20-1410T to 20-1450T

16 MM 20-1610T to 20-1650T

18 MM 20-1810T to 20-1850T

►LENGTH IN MM 10, 15, 20, 25, 30, 35, 40, 45, 50



EXPANDABLE CAGE

►DIA Ø ►REF. NO. TT

12 MM 20-1225T to 20-1277T

14 MM 20-1425T to 20-1477T

16 MM 20-1625T to 20-1677T

18 MM 20-1825T to 20-1877T

►LENGTH IN MM 20 X 25, 25 X 35, 35 X 55, 55 X 77



SPINAL IMPLANTS

TLIF CAGE



►DEGREE°	►REF. NO. TT
4°	20-0407T to 20-0411T
8°	20-0807T to 20-0811T
11°	20-1107T to 20-1111T
►LENGTH IN MM	7 MM X 25 MM, 8 MM X 25 MM
	9 MM X 25 MM, 10 MM X 25 MM
	11 MM X 25 MM

PLIF CAGE



►DEGREE°	►REF. NO. TT
0°	20-0007T to 20-0014T
4°	20-0407T to 20-0414T
►LENGTH IN DIA	7Ø, 8Ø, 9Ø, 10Ø, 11Ø, 12Ø, 13Ø, 14Ø

ROD TO ROD CLAMP

ROD TO ROD CLAMP - AXIAL



ROD TO ROD CLAMP - PARALLEL



►PRODUCT	►REF. NO. TT
ROD TO ROD CLAMP - AXIAL	20-0001T
ROD TO ROD CLAMP - PARALLEL	20-0002T

CERVICAL PLATE STABILISATION PIN



►REF. NO. TT
20-0001T

SPINAL IMPLANTS

TRANSVERS CONNECTOR



►DIA Ø ►REF. NO. TT
5.5 MM 20-5500T

CERVICAL SPINE OFFSET



►REF. NO. TT
20-0000T

CERVICAL ROD CONNECTOR



►REF. NO. TT
20-0000T

HALF PIN



►DIA Ø ►REF. NO. TT
5.0 MM 20-50150T
►LENGTH IN MM
150 MM



Thanks For Going Throug Catalog

COMPANY PROFILE:

Company Name:

Sinus instruments .co

Establishment

1986

Business:

Manufacturers-cum-Exporters of
Surgical, Dental, Beauty Instruments
& all sort of Scissors

Telephones:

0092523301054

Cell:

+923049315806

Sinusinstruments@gmail.com